

# Operating conduino through a MATLAB script

Novak Lab, Newcastle University

September 6, 2025

Conduino is a high-resolution multichannel water conductivity sensor. It uses micro-USB connectors for fast sensing of bipolar or tetrapolar liquid conductivity with minimal spatial granularity. The subsequent sections discuss a complete procedure for collecting data from the conduino board through a MATLAB script. Detailed information on conduino can be found [here](#).

## 1 Getting the codes

Compiling and uploading [Firmware Conduino Matlab.ino](#) code is required to upload the arduino sketch to the conduino board. After uploading the arduino sketch to the board, the data can be collected from the controller board using the MATLAB script [ConduinoMultichannel.m](#).




Figure 1: Arduino IDE interface

## 2 Uploading Arduino sketch on conduino board

An open source software [Arduino IDE](#) should be installed to upload arduino sketch ([Firmware Conduino Matlab.ino](#)) on the conduino board. An administration log-in is required for the software installation, compilation, and uploading of the code. After successful installation of the software, conduino board should be connected to the PC. A USB port powers the board, so there is no need to be a separate connection for the power supply. After connecting the board, open the ([Firmware Conduino Matlab.ino](#) code on Arduino IDE software. Figure 1 shows the page that will appear when you open the ([Firmware Conduino Matlab.ino](#) code. Select the *Arduino Uno* board, which will automatically appear after clicking the select board option. The serial port number will also appear with the *Arduino Uno* board (Figure 1). To compile the code, click the compile icon shown in Figure 1. After successful compilation, the code can be uploaded to the board by clicking the upload icon shown in Figure 1.




Figure 2: Step 1

## 3 Data collection through MATLAB script

To operate conduino through MATLAB, an add-on package *MATLAB Support Package for Arduino Hardware* should be installed from the *Add-Ons* tab of the *Environment* section from the MATLAB toolstrip (Figure 2). An administration log-in is required to complete the installation, and conduino board should be connected to the PC via USB port. After installing the package, a prompt will appear to proceed with the configuration. The configuration can be completed by following the instructions indicated in Figures 3, and 4. In Figure 4, the board and the port should be selected correctly; the correct board and port option will be prompted during selection. The libraries *I2C*, *SPI*, and *Servo* should be selected, and libraries are uploaded to the server by clicking the *program* tab. The *Test connection* tab should be clicked to check if all the steps

are completed correctly (Figure 5). After a successful test, the configuration can be completed by clicking *Next*.




Figure 3: Step 2




Figure 4: Step 3

After the successful completion of configuration, and with the board connected to the PC, open the [ConduinoMultichannel.m](#) script on MATLAB. Remember to change the port number in the code. The correct port number can be found by checking the COM port in the Arduino IDE interface shown in Figure 1. The data can be collected by running the script.




Figure 5: Step 4

## 4 Protocol for plotting background stratification

To obtain background stratification [Stratification.m](#) should be run. The code computes the density variation along the height moved by the rig. The steps to be followed to collect data from the probe are;

Step 1: Place the rig above the top layer of the water with its tip touching the water surface.

Step 2: Connect the Conduino Board.

Step 3: Select the appropriate port ( refer to section 3).

Step 4: Run the code [Stratification.m](#). In this step threshold value for the data to be filtered out should be determined from *Figure 1.fig*. The data below the data point corresponding to the circled portion in Fig. 6 should be filtered out.

Step 5: Run the code [Stratification.m](#) again after including the threshold data value in the code.

Step 6: Enter the number of samples ( e.g. 300) and start the motor by pressing the green button on the stepper motor firmware.

Step 7: Save the *Stratification.fig* file generated after collecting data from probes in the folder containing the [Stratification.m](#) and [Processed\\_Image\\_Tangent.m](#). Collect the data stored in array `ch1[ ]`, `C[ ]`, `rho[ ]`, `d[ ]`, and `h[ ]` and paste it in the Excel file (*Probe Reading.xls*) shown in Fig. 7.

Step 8: Run the code [Processed\\_Image\\_Tangent.m](#) to calculate the pycnocline thickness. The code generates a tangent on *Stratification.fig* and draws a horizontal line through the tangent point and calculates the distance between it.

Step 9: Save the *Sratification.fig*, *Raw Data.fig*, *Smooth\_Stratified.fig* (if needed) and *Probe Reading.xls* in the folder *DDMMYY* that also contains digiflow data.




Figure 6: Threshold Data

**Note:** In step 1, placing the probe tip just at the surface of the top layer gives the start signal to filter out the data required from the raw data collected by the probe. As data only needed to be collected during the forward cycle of the rig movement, the initial rig placement also gives stop signal for collecting the required data. In step 5, running the code [Stratification.m](#) will collect data by the rig during the forward cycle only. The detailed explanation about the code is explained in the section 6. The code will generate three figures *Figure 1.fig*, *Figure 2.fig* and *Figure 3.fig*. Save the *Figure 1.fig* as *Raw Data.fig* (refer to step 9), *Figure 2.fig* as *Stratification.fig* (refer to step 7 and step 9), and *Figure 3.fig* as *Smooth\_Stratified.fig* (refer to step 9). It is to be further noted that saving *Figure 3.fig* is only required when there is visible noise in the *Figure 2.fig*. Furthermore, in step 7, the *Stratification.fig* is stored in the folder with the code [Processed\\_Image\\_Tangent.m](#) for image processing. The generated image will give information about the pycnocline thickness

Probe Readings	
Date:	
Number of Samples:	
Time interval between Samples (ms):	
Observation 1	
Raw Data_Channel 1 (Ch1[]):	
Channel 1_Trimmed Data (C[]):	
Channel 1_Trimmed Data (rho[]):	
Distance Travelled by Rig (d[]):	
Distance Travelled by Rig converted to height (h[]):	
Observation 2	
Raw Data_Channel 1 (Ch1[]):	
Channel 1_Trimmed Data (C[]):	
<div> <div>&lt;</div> <div>&gt;</div> <div>Main_Tank</div> <div>Behind_Gate</div> <div>+</div> </div>	

Figure 7: Probe Reading.xls

## 5 Code for changing the rig travel pattern

The distance moved by the rig can be modified by uploading code [Stepper\\_Motor.ino](#) and [Stepper\\_motor\\_5mm\\_step.ino](#) on the firmware that moves the stepper motor. To upload the code on the stepper motor firmware, connect the firmware to the computer via USB connection. It is to be noted that the rig should not be connected to other power source during the upload. The firmware derives power from the computer via USB connection; if the rig is connected to another power source, it will damage the firmware, as it will receive power from two sources. Open the code [Stepper\\_Motor.ino](#) on [Arduino IDE](#) interface to modify the rig travel distance. The total travel distance can be modified by changing line 27 of the code [Stepper\\_Motor.ino](#) highlighted in Fig. 8a. At present, it is configured to travel a distance of 240 mm (Fig. 8a). The rig movement can also be modified to travel in small steps by uploading the code on the stepper motor firmware. The step distance can be modified by modifying line 99 of the code as highlighted in Fig. 8b. The code now will move the rig in steps of 5 mm (Fig. 8b). It is to be noted that the firmware can be configured with only one code at a time.

```
17 // Start and Stop LED controls
18 const int startLED = 2; //Start LED pin
19 const int stopLED = 3; //Stop LED pin
20 int startLEDState = LOW; //Start LED State
21 int stopLEDState = LOW; //Stop LED State F
22 int startButton = 0; //Start button val
23 int stopButton = 0; //Stop button valu
24 int limit = 0; //Initial limit co
25 unsigned long int stepsPerRevolution = 24550
26
27 int DIST = 240; //Distance in mm {
28
```

(a)

```
91
92 void MoveStepper() {
93     unsigned long currentMillis = millis();
94     startButton = analogRead(startBTN);
95     // Serial.println(startButton);
96     if (startButton <= 500 && limit >=250){
97         digitalWrite(sampletrig, HIGH);
98         delay(50);
99         stepper.runRelative(-5);
100         delay(50);
101         //digitalWrite(sampletrig, LOW);
102     }
103     delay(1);
104
105 }
```

(b)

Figure 8: Stepper Motor Configuration

## 6 Matlab code for computation of background stratification

The code [Stratification.m](#) computes the probe data converted to the density reading along the vertical height of the tank. The theory behind the calculation of distance travelled by the rig is discussed in the section 8. In the present code, the rig travels 24 cms of total distance in the forward cycle, where the initial 1.2 cms it travels with acceleration given by equation 1. However, the data collected by probe during this period is filtered out, thus distance calculation during this period is not included in the code. The different portions of the code are explained in blocks as follows.

BLOCK 1 (Fig. 9) : The code takes input like port number and system type ( mac or pc) from the user. It also defines the acceleration and velocity values determined by checking the time ti from the code (refer to section 8). In the block, different arrays are defined to store the time information from the probe in terms of filtered data. The array ch1[] stores all the raw data from the probe, while cht[] only stores data above threshold values. The threshold data is determined during step 4 of the section 4. Furthermore, chf[] stores the first half of the filtered out data from cht[]. The vartime [] stores the corresponding time of the cht[]. while recordtime [] will start the time counter from zero by taking readings from vartime []. th [] stores the first half of the trimmed data from recordtime[] that also considers the initial 1.2 cms of rig movement in account,

while `thf []` will start the time counter from zero by taking readings from `th []`. The array `d []` and `t []` is calculated according the information given in equations 1, 2 and 3. The `C[]` stores the final trimmed value of conductivity readings that is required to plot the stratification. `rho []` stores the converted density values.

```

1  %stratification.m
2  clear
3  close all
4  delete(instrfindall); % Use this if the ports get stuck
5
6  %%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%% Settings %%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%
7  macORpc = 'pc'; % 'mac' or 'pc'
8  portNumber = 10; % COM port number on pc, port index on mac
9  a1 = 10.41; %Motor acceleration and deceleration value
10 v1 = 5; %Motor Constant velocity
11 vartime = []; %array to store the time values corresponding to the conductivity read
12 recordtime = []; %array that will start the time counter based on vartime
13 %d1 = []; %array to store the distance travelled by rig during acceleration
14 d2 = []; %array to store the distance travelled by rig with constant velocity
15 d3 = []; %array to store the distance travelled by rig during deceleration
16 d = []; %array to store the total distance travelled by rig
17 cht = []; %array to store the conductivity readings above threshold value
18 chf = []; %array to store the the first half of the conductivity readings above thre
19 th = []; %array to store the the time corresponding to the first half of the conduct
20 thf = []; %array that will start time counter based on th[]
21 %t1 = [];
22 t2 = []; %array to store the thf[] values when rig travels with constant velocityc
23 t3 = []; %array to store the thf[] when rig travels during deceleration
24 t = []; %array to store the t2[] and t3[] values
25 h = []; %array to store to convert the distance values from d[] to height starting
26 %htf = [];
27 C = []; %array that will store the values of the conductivity readings based on dist
28 rho = []; %array that will store the vdensity readings corresponding to the conduct
29 %%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%%
30

```

Figure 9: BLOCK 1

BLOCK 2 (Fig. 10): In Block 2 of the code, the code reads raw data from the serial port and stores the values in the arrays of the connected channels. In most cases, channel 1 is connected, and raw data is stored in `ch[]`.

BLOCK 3 (Fig. 11): In Block 3 of the code, the data filtering is done, where the user gives the value of the threshold data below which all the data should be filtered out. The portion of the code also computes the first half of the probe data.

BLOCK 4 (Fig. 12): In Block 4 the distance computation are done. It is done by running two loops. each for the calculation of `d2[]` and `d3[]` respectively ( equations 2, 3), The loop starts with `t1` and ends at  $9*t1$  for `d2[]`, while for `d3[]` it ends with  $10*t1$ .

BLOCK 5 (Fig. 13): In Block 5 the conductivity reading is stored in array `C[]` for plotting against `h[]`. The `h[]` computation is also done during this block of the code. The conversion of conductivity readings from `C[]` to `rho[]` is also done in this block of the code. The conversion is done according to the equation 9.

BLOCK 6 (Fig. 14): Sometimes the obtained conductivity data has noise; in order to smooth the curve of the plot, data smoothing can be done if required. This is done in block 6 of the code. In this block, code to detect the point for drawing tangent through the curve is also added.

BLOCK 7 (Fig. 15): In this block, code to generate three figures discussed in step 8 of the section 4 is mentioned. It is to be noted a separate code [Processed\\_Image\\_Tangent.m](#) for drawing tangent and getting information about pycnocline thickness is also available. The tangent and thickness calculation by this code will appear as Fig. 16.

```

%assigns the object s to serial port
switch macORpc
    case 'pc'
        s1 = serial(['COM' num2str(portNumber)]); % for PC
    otherwise
        s = seriallist;
        disp(['Using port ' s(portNumber) ])
        s1 = serial(s(portNumber)); % for mac
end

set(s1, 'InputBufferSize', 128); %number of bytes in inout buffer
set(s1, 'BaudRate', 115200);
set(s1, 'Parity', 'none');
fopen(s1);
clc
prompt = 'How many samples? (dt ~40 ms typically) ';
nSamples=input(prompt,'s');
fprintf(s1,'%c',nSamples);
nSamples = str2double(nSamples);
ch1=zeros(0,nSamples);
ch2=zeros(0,nSamples);
ch3=zeros(0,nSamples);
ch4=zeros(0,nSamples);

tic
fwrite(s1, '%f');
for nSample=1:nSamples
    ch1(nSample)= fscanff(s1, '%f');
    ch2(nSample)= fscanff(s1, '%f');
    ch3(nSample)= fscanff(s1, '%f');
    ch4(nSample)= fscanff(s1, '%f');

```

BLOCK 2

Figure 10: BLOCK 2

```

for i = 1:length(time) %loop for storing conductivity readings above threshold value and the corresponding time
    if ch1(i) > 3.1 %threshold need to be confirmed from Figure1
        vartime = [vartime, time(i)];%it stores the successive values of vartime
        recordtime = [vartime - vartime(1)];%it stores the successive values of recordtime
        cht = [cht, ch1(i)];%it stores the successive values of conductivity readings
    end
end

halflength = floor(length(cht)/2); %code to store the first half of the conductivity readings and the corresponding vartime
chf = cht(1:halflength);
th = vartime(1:halflength);
for j = 1:length(thf)
    %thf = [thf, th(j)-th(1) + 0.516];%code to store the successive values of thf
    thf = [thf, th(j)-th(1)];%code to store the successive values of thf
end
ti = thf(end)/10;%the rig moves with deceleration for time ti.

```

BLOCK 3

Figure 11: BLOCK 3

## 7 Appendix

### 8 Rig travel distance calculation

The rig driven by a stepper motor travels with an initial 5 % of the distance with constant acceleration ( $a$ ), next 90 % with constant velocity, and the last 5 % with constant deceleration. Thus, to calculate the distance ( $d$ ) in time ( $t$ ), let us consider the rig travels a distance of  $d_1$  in time  $t_1$  with acceleration,  $d_2$  in time  $t_2$  with constant velocity ( $v$ ),  $d_3$  in time  $t_3$  with constant deceleration ( $a$ ). Thus we will get,



```

100 % the constant velocity
101 starttime = thf(1); %loop starts from initial value of thf
102 endtime = 9*ti; %loop ends at time when rig travelled with constant velocity
103 %maxIter = 13;
104 i = 1; % starting loop counter
105 while thf(i) >= starttime && thf(i) < endtime %loop will run during the constant velocity run of rig
106     d2 = [d2, 1.2 + (thf(i)-thf(1))*v1]; %distance calculation based on constant velocity, the threshold conductivity reading start
107     t2 = [t2, thf(i)];
108     i = i + 1;
109 end
110 % Portion of the code to calculate the distance travelled by rig during
111 % the deceleration of the motor
112
113 starttime = 9*ti; %loop starts from value of thf when rig starts movement with constant deceleration
114 endtime = 10*ti; %loop ends when rig completes first half of the travel.
115 %maxIter = 13;
116 i = i(end); %loop starts counter right after the previous loop counter ends
117 while thf(i) >= starttime && thf(i) < endtime %loop during the deceleration run of rig
118     d3 = [d3, d2(end) + (thf(i)-t2(end))*v1 - (thf(i)-t2(end))*(thf(i)-t2(end))*0.5*a1];
119     t3 = [t3, thf(i)]; % distance calculation during deceleration run of rig
120     i = i + 1;
121 end
122 d = [d2 d3];

```

Figure 12: BLOCK 4

```

123 starttime = thf(1); %loop to filter the conductivity reading values based on distance calculation
124 endtime = t3(end); %the end time should be taken from the t3[] in order to remove any error arising due to data trimming
125 %maxIter = 13;
126 j = 1;
127 while thf(j) >= starttime && thf(j) <= endtime
128     C = [C, chf(j)];
129     h = [h, (41.2 - d(j))]; %it will convert the rig travel distance to the height moved by rig in the tank starting from 40 cm
130     j = j + 1;
131 end
132 t = [t2 t3]; % array to store time values from thf[] with size equal to C[];
133 starttime = thf(1); %loop to convert the conductivity reading to the corresponding density values
134 endtime = t3(end); %the end time should be taken from the t3[] in order to remove any error arising due to data trimming
135 %maxIter = 13;
136 k = 1;
137 while thf(k) >= starttime && thf(k) <= endtime
138     rho = [rho, (1025 + ((C(k) - C(1))/(C(end) - C(1)))*25)]; %linear interpolation of data for converting conductivity readings
139     k = k + 1;
140 end

```

Figure 13: BLOCK 5

```

145 rho_smooth = smoothdata(rho, 'sgolay', 15);
146 dh = gradient(h, rho_smooth);
147 d2h = gradient(dh, rho_smooth);
148 % it will detect the index of the point where tangent to be drawn
149 idx_inflect = find(diff(sign(d2h))); % it will detect the index of the point where tangent to be drawn
150
151 %dh = gradient(h, rho);
152 %d2h = gradient(dh, rho);
153 % Choose the index of the point where you want the tangent
154 %idx_inflect = find(diff(sign(d2h))); % Change this index as needed (between 2 and length(x)-1)
155
156
157
158 plot(ch1); % Plots the total number of samples (figure1)
159 hold on;
160 plot(ch2);
161 plot(ch3);
162 plot(ch4);
163 grid on
164 xlabel('Sample')
165 ylabel('Reading')
166

```

Figure 14: BLOCK 6

$$d_1 = 0.5at_1^2 \rightarrow 0.05d = 0.5at_1^2 \rightarrow a = \frac{0.1d}{t_1^2} \text{ and } v = at_1 \rightarrow v = \frac{0.1d}{t_1} \quad (1)$$

$$d_2 = vt_2 \rightarrow 0.9d = \frac{0.1d}{t_1}t_2 \rightarrow 9t_1 = t_2 \quad (2)$$

$$d_3 = vt_3 - 0.5at_3^2 \rightarrow 0.05d = \frac{0.1d}{t_1}t_3 - 0.5\frac{0.1d}{t_1^2}t_3^2 \rightarrow 1 - 2\frac{t_3}{t_1} + \frac{t_3^2}{t_1^2} = 0 \rightarrow t_3 = t_1 \quad (3)$$

Thus, it can be found out that the rig travels a total distance  $d$  in total time  $11t_1$ . The total rig travel distance can be set by using the codes given in section 5. The time  $t_1$  value can be found by dividing the  $\text{thf}[\text{end}]$  ([Stratification.m](#)) by 11. In the code [Stratification.m](#), the  $t_1$  is calculated as  $t_i$ , which is calculated by dividing  $\text{thf}[\text{end}]$  by 10. As in the code [Stratification.m](#) the initial acceleration period is not considered this is why the  $\text{thf}[\text{end}]$  value is divided by 10.

```

198 plot(rho_smooth,h);%plots the variation of smoothed density with the height
199 hold on;
200 %plot(time,ch2);
201 %plot(time,ch3);
202 %plot(time,ch4);
203 grid on
204 for n = 1:length(idx_inflect)%plots the tangent in the same figure
205     i0 = idx_inflect(n);
206     rho0 = rho_smooth(i0);
207     h0 = h(i0);
208     slope = dh(i0);
209
210     tan_rho = linspace(rho0 - 0.5, rho0 + 0.5, 100);
211     tan_h = h0 + slope*(tan_rho - rho0);
212     %plot(rho0, h0);
213     %plot(rho0, h0, 'ro', 'MarkerSize', 8, 'LineWidth', 2);
214     %plot(tan_rho, tan_h, 'r--', 'LineWidth', 1.5);
215     plot(tan_rho, tan_h);
216     %plot([min(rho) max(rho)], [h0 h0]); %plots and horizontal line through
217     %the tangent point.
218
219 end
220 xlabel('Density (kg/m^3)')
221 ylabel('Height')
222 %legend('ch1', 'ch2', 'ch3', 'ch4')
223 pp=axis;
224 axis([1000 1060 10 50])
225

```

BLOCK 7

Figure 15: BLOCK 7




Figure 16: Code for Image Processing

## 9 Linear interpolation formula

The density information corresponding (values in the array  $\rho_0$ ) to the conductivity readings (values in the array  $C$ ) can be obtained from linear interpolation. The corresponding density values are obtained by the given relation,

$$\rho_i = \rho_0 + \frac{(\rho_n - \rho_0)}{(C_n - C_0)}(C_i - C_0) \quad (4)$$

where  $\rho_0$ ,  $\rho_i$ , and  $\rho_n$  are the density of the upper layer, the  $i^{th}$  layer and the lower

layer.  $C_0$ ,  $C_i$  and  $C_n$  are the first element,  $i^{th}$  element, and last element of  $C[ ]$ .